

Week 26

Monday 26/6	3
Verification of Cas9-STE2-STE3 vector	3
Purification	3
Restriction analysis	4
Yeast transformation STE2gRNA and STE3gRNA, and Raphaels Cas9 plasmid with gRNA	6
Preparation of medium and O/N-culture	6
PCR amplification of Cre-Cas9 constructs	6
Preparations and PCR	6
Gel electrophoresis	7
PCR purification	7
Water evaporation of GPCRs	7
Evaporation for increasing the concentration of the GPCRs	7
Tuesday 27/6	8
Verification of Cas9-STE2-STE3 vector	8
Sequencing	8
Yeast transformation STE2gRNA and STE3gRNA, and Raphaels Cas9 plasmid with gRNA	9
Transformation	9
Gel electrophoresis and purification of PCR fragments Cre-Cas9	10
Gel electrophoresis	10
PCR amplification of Cre-Cas9 constructs	11
Preparation and PCR	11
Gel electrophoresis	11
Gel extraction	11
PCR amplification of Cre-gRNA constructs	12
Preparation and PCR	12
Gel electrophoresis	12
PCR amplification of Cre-gRNA constructs and Cre-Cas9_5 construct	12
Preparation and PCR	12
Gel electrophoresis	13
PCR purification of gRNA fragments 1 and 5.	13
Fusion PCR of CreCas9s construct	13
PCR	13
Wednesday 28/6	15
Fusion PCR of CreCas9s construct	15
PCR	15
Date: 26/6	17
What was done	17
	1

PCR amplification of Cre-gRNA constructs	17
Preparation and PCR	17
Gel electrophoresis	18
PCR Purification of gRNA_2	18
Date: 26/6	18
What was done	18
PCR amplification of Cre-gRNA constructs	18
Preparation and PCR	18
Gel electrophoresis	18
Thursday 29/6	20
PCR amplification of Cre-gRNA constructs 1,3 and 4.5 and Cre-Cas9 construct 5	20
Preparations and PCR	20
Gel electrophoresis	20
Date: 29/6	21
What was done	21
PCR amplification of Cre-gRNA constructs 1,3 and 4.5	21
Preparations and PCR	21
Gel electrophoresis	21
Date: 29/6	21
What was done	21
Fusion PCR of CreCas9s construct	21
Gel electrophoresis	21
Date: 29/6	22
What was done	22
Gibson assembly for Cas9 with p416 vector	22
Date: 29/6	23
What was done	23
PCR amplification of repair fragments O1, O2, R1 and R2	23
Preparations and PCR	23
Yeast transformation of STE2gRNA and STE3gRNA, and Raphaels Cas9 plasmid with gRNA	24
O/N-culture	24
Friday 30/6	25
Gel Extraction of gRNA fragment 3	25
Date: 30/6	25
What was done	25
Gel electrophoresis of PCR products O1, O2, R1 and R2	25
Date: 30/6	25
What was done	25
PCR amplification of repair fragments O1, O2, R1 and R2 as well as gRNA fragment 4.5	25

Gel electrophoresis on GPCRs and gRNA 4.5	26
Yeast transformation of STE2gRNA and STE3gRNA, and Raphaels Cas9 plasmid with gRNA	26
Transformation	26
Date: 30/6	27
What was done	27
Touchdown PCR of gRNA 4.5	27
Preparations and PCR	27


** Template → More detail written how to do stuff (for example the first time we do it)

Monday 26/6

Date: 26/6	Performed by: Alex & Angelo																						
What was done	Observations/comments																						
Verification of Cas9-STE2-STE3 vector																							
Purification																							
<p>The 6 tubes from each strain (STE2 and STE3) was purified according to the purification protocol (plasmid miniprep kit). The following concentration was observed after the purification.</p> <p>The colonies from STE2 was marked with X.2, where X is the colony, and STE3 was marked with X.3.</p> <p>Those which are marked in red are those which were sent in for sequencing. Based on the concentration and which were clearest to see on the gel. See next day (27/6) for more information about sequencing.</p>	<p>Concentration:</p> <table> <tr> <th>Tube</th><th>conc. [ng/μL]</th></tr> <tr> <td>1.2</td><td>208,5</td></tr> <tr> <td>2.2</td><td>295,1</td></tr> <tr> <td>2.3</td><td>272,3</td></tr> <tr> <td>4.2</td><td>357,7</td></tr> <tr> <td>5.2</td><td>399,7</td></tr> <tr> <td>6.2</td><td>273,9</td></tr> </table> <table> <tr> <th>Tube</th><th>conc. [ng/μL]</th></tr> <tr> <td>1.3</td><td>369,5</td></tr> <tr> <td>2.3</td><td>290,1</td></tr> <tr> <td>3.3</td><td>335,8</td></tr> </table>	Tube	conc. [ng/μL]	1.2	208,5	2.2	295,1	2.3	272,3	4.2	357,7	5.2	399,7	6.2	273,9	Tube	conc. [ng/μL]	1.3	369,5	2.3	290,1	3.3	335,8
Tube	conc. [ng/μL]																						
1.2	208,5																						
2.2	295,1																						
2.3	272,3																						
4.2	357,7																						
5.2	399,7																						
6.2	273,9																						
Tube	conc. [ng/μL]																						
1.3	369,5																						
2.3	290,1																						
3.3	335,8																						

	<table> <tr> <td>4.3</td><td>483,9</td></tr> <tr> <td>5.3</td><td>360,9</td></tr> <tr> <td>6.3</td><td>365,6</td></tr> </table>	4.3	483,9	5.3	360,9	6.3	365,6		
4.3	483,9								
5.3	360,9								
6.3	365,6								
<h3>Restriction analysis</h3>									
<p>The digestions were done according to the FastDigest protocol. A master mix was made (x15)</p> <p>210 µL Water 30 µL FastDigest buffer 15 µL NheI 15 µL BamHI</p> <p>18 µL of the master mix was mixed with 2 µL DNA (the samples).</p>	<p>2 µL DNA was calculated from 4.3. The max amount of DNA should be < 1 µg. So the concentration for 4.3 is 483,9 ng/µL → 1000 ng / 483,9 ng/µL = 2 µL. The same volume was used for all digestions.</p>								
<p>2 µL was loaded in every well (The bands became messy) 0,5 µL of ladder was added</p> <p>1 µL was loaded in the second try, and the bands became much clearer.</p> <p>In the second try a misload appeared for colony 5, and was loaded twice.</p> <p>The samples sent for sequencing were 4.2 and 6.2 for STE2 and for STE3 5.3 and 6.3. See text below!</p>	<p>Expected results for the restriction analysis can be seen below</p> <p>The expected result is based on the program Benchling's approximation for digestion. The control (backbone plasmid Cas9-vector) is supposed to give two bands. If the insert had been successfully inserted in the backbone plasmid three bands should be visible.</p> <div> <table> <tr> <td>Ladder</td><td>GeneRuler 1 kb Plus</td></tr> <tr> <td>1</td><td>Cas9 with gRNA STE2 vector - BamHI NheI</td></tr> <tr> <td>2</td><td>Cas9 vector - BamHI NheI</td></tr> <tr> <td>3</td><td>Cas9 with gRNA STE3 vector - BamHI NheI</td></tr> </table> </div>	Ladder	GeneRuler 1 kb Plus	1	Cas9 with gRNA STE2 vector - BamHI NheI	2	Cas9 vector - BamHI NheI	3	Cas9 with gRNA STE3 vector - BamHI NheI
Ladder	GeneRuler 1 kb Plus								
1	Cas9 with gRNA STE2 vector - BamHI NheI								
2	Cas9 vector - BamHI NheI								
3	Cas9 with gRNA STE3 vector - BamHI NheI								

The results from restriction analysis is shown to the left. The Lad (Ladder) used was GeneRuler 1 kb. Well 1 is the control (backbone plasmid cut with RE) and well 2-7 is our samples (the different colonies). All of the colonies for STE2 showed results according to the prediction, and for STE3 all except one (colony 4 seemed a bit off) produced expected results. Based on this and the concentration shown in Verification (27/6), purification step above and those which looked best on the gel sent for sequencing.



Date: 26/6	Performed by: Alex & Angelo									
What was done	Observations/comments									
Yeast transformation STE2gRNA and STE3gRNA, and Raphaels Cas9 plasmid with gRNA										
Preparation of medium and O/N-culture										
Medium that's need to be autoclaved: Sorbitol 1M Preparation of 50 mL: 50 mL MQ + 9,1 g Sorbitol LiAc 1M Preparation of 50 mL: 50 mL MQ + 5,101 g LiAc PEG 50% w/v Preparation of 50 mL: 50 mL MQ + 25 g PEG O/N-culture: 2 yeast culture was inoculated in 5 mL YPD along with a negative control (falcon 50 mL tubes). One colony was picked from each strain. The O/N-culture were put in a shake room 30°C overnight.										
<table><tr><th>Yeast strain</th><th>deletion</th><th>Mating type</th></tr><tr><td>11C</td><td>$\Delta HIS, \Delta URA$</td><td>a</td></tr><tr><td>61A</td><td>$\Delta HIS, \Delta URA, \Delta LEU$</td><td>$\alpha$</td></tr></table>	Yeast strain	deletion	Mating type	11C	$\Delta HIS, \Delta URA$	a	61A	$\Delta HIS, \Delta URA, \Delta LEU$	α	OBS! Important to remember which is a and α !!
Yeast strain	deletion	Mating type								
11C	$\Delta HIS, \Delta URA$	a								
61A	$\Delta HIS, \Delta URA, \Delta LEU$	α								

Date: 26/6	Performed by: Group 1
What was done	Observations/comments
PCR amplification of Cre-Cas9 constructs	
Preparations and PCR	
<p>IDT blocks and primers for the Cre-Cas9 constructs and Cre-gRNA were diluted 1:10 in their original tubes.</p> <p>The constructs and primers for the Cre-Cas9 was diluted again before the PCR.</p> <p>The PCR was done according to the standard Phusion protocol with an annealing temperature at 60°. The constructs used were the Cas9 constructs 1, 2, 3, 4, 4.5, 5 and "middle". Middle was the long Cas9 part and used a plasmid as template.</p>	<p>Note: Some of the constructs had a original amount of 500 ng and therefore the dilution became 5 ng/μl before this was discovered. This is noted on the tubes!</p> <p>The second dilutions were all made correctly so that the concentrations was 1 ng/μl.</p>
Gel electrophoresis	
<p>A gel electrophoresis was performed.</p> <p>The constructs loaded were in the order: 1, 2, 3, 4, 4.5, 5 and "middle".</p>	<p>1 - wrong - remake 2 - ok - PCR purification 3 - wrong - remake 4 - ok - PCR purification 4.5 - ok - PCR purification 5 - wrong - remake middle - ok?? - digested - PCR purification</p>
PCR purification	
<p>PCR purification was performed with the GeneJET purification kit.</p> <p>Samples were the gRNA constructs 2, 4, 4.5 & middle.</p> <p>Concentrations were then checked with nano drop.</p>	<p>Concentrations were:</p> <p>2: 0,3 ng/μl. - remake 4: 4,8 ng/μl. - remake 4.5: 21,7 ng/μl - ok! middle: 21 ng/μl - ok!</p>

Date: 26/6	Performed by: Group 1
What was done	Observations/comments
Water evaporation of GPCRs	
Evaporation for increasing the concentration of the GPCRs	
<p>Two samples each of the constructs O1 (Olf1258 part 1) O2, R1 (Rat17 part 1) and R2 were mixed and evaporated together in the evaporation machine.</p> <p>After the evaporation the DNA from O2, R1 and R2 were resuspended in 20 µl water. O2 wasn't fully evaporated.</p> <p>The concentrations were measured with nano drop.</p>	<p>The volumes of the samples after the resuspension and after nanodrop was:</p> <p>O1: 16 µl O2: 19 µl R1: 19 µl R2: 19 µl</p> <p>The concentrations were:</p> <p>O1: 126,3 ng/µl O2: 152,9 ng/µl R1: 199,6 ng/µl R2: 80,8 ng/µl</p> <p>Total milligram:</p> <p>O1: 2,02 O2: 2,91 R1: 3,79 R2: 1,54</p>

Tuesday 27/6

Date: 27/6	Performed by: Alex & Angelo
What was done	Observations/comments
Verification of Cas9-STE2-STE3 vector	
Sequencing	
<p>The total volume, V_{Total} should be 15 µL, and the concentration, c_{DNA} of DNA 50 ng/µL.</p> <p>$c_1 \cdot V_1 = c_2 \cdot V_2$</p> <p>$c_{DNA} \times V_{Total} (15 \mu L) = c_{sample} \times V_{Need to be added}$</p>	$\frac{50 \frac{ng}{\mu L} \times 15 \mu L}{c} = V_{Need to be added}$ <p>The volume MQ needed was calculated by</p> <p>$V_{Total} - V_{Need to be added}$</p>

With this equation the volume needed of our purified plasmids was calculated as following:

Sample	Conc. [ng/ μ L]	Volume DNA	Volume MQ
STE2 colony 4 (4.2)	357,7	2,7	12,3
STE2 colony 6 (6.2)	273,9	2,1	12,9
STE3 colony 5 (5.3)	360,9	2,1	12,9
STE3 colony 6 (6.3)	365,6	2,1	12,9

4 tubes were prepared. 4 stickers with following labelled was taken and match up

AAA0085 073 4.2

AAA0085 074 4.6

AAA0085 075 5.3

AAA0085 076 6.3

AAA0085 077 Primer

The primers are designed on the computer later on when you will sent in the order. A supervisor will help you with this step.

Fast way to calculate is always to take 750 divided by the concentration of the plasmid or construct yo will sent in ($15 * 50 = 750$).

Keep in mind that the following with highest concentration were chosen in the transformation, marked in red.

You will find the stickers in a folder in the large lab on the other side of the MQ-water tap, across the refrigerator. Were the risk declaration will be handed in. Write down which stickers you take (for example AAA0085 73-77) and sign.

Date: 27/6	Performed by: Alex & Angelo
What was done	Observations/comments
<p align="center">Yeast transformation STE2gRNA and STE3gRNA, and Raphaels Cas9 plasmid with gRNA</p>	
<p align="center">Transformation</p>	
<p>The two O/N-cultures was diluted 10x (45 mL YPD was added to the 50 mL falcon tubes and then transferred to a shake flask, one for each strain). A OD measurement were made at 8:15 for a 100x dilution. The shake flasks was put in shaker room for 2 h 20 min. →</p>	<p>Following OD was measured</p> <p align="center">100x 10x</p> <p>11C ~ 0,03 → 0,3</p> <p>61A ~ 0,02 → 0,2</p> <p>→ This may differ, see protocol for how long. It depends on OD.</p>

<p>After the 2 h 20 min the LiAc protocol was followed for the whole transformation.</p> <p>After the heat shock the tubes were centrifuge and the supernatant was discarded. A recovery step was added: 500 mL Sorbitol and 500 mL YPD was added. The tubes were put in incubation at 30°C for 1 h. Then the protocol was followed again.</p>	<p>NOTE! Some steps were made different. When the PEG, LiAC and SS-DNA was added the order was not the same as the protocol said. And the incubation for 30°C was skipped. Some other rfc was used than what the protocol says and different times.</p> <p>The steps that we skipped might be bad (have not seen the results yet) but follow the protocol instead for skipping steps and so on to be sure. The recovery step was good due!</p>
--	--

Date: 27/6	Performed by: Group 1
What was done	Observations/comments
<h2 style="text-align: center;">Gel electrophoresis and purification of PCR fragments Cre-Cas9</h2>	
<h3 style="text-align: center;">Gel electrophoresis</h3>	
<p>The fragments were loaded in the following order: 1, 2, 3, 4, 5. The gel was run for 35 min at 75 V, 160 mA.</p> <p>A picture was taken on the gel with UV-light.</p> <hr/> <p>GEL EXTRACTION First the sample was prepared for gel extraction. 8 µl loading dye was added to 48 µl PCR-product. It was then loaded on the gel.</p> <p>Gel extraction was done on Cre-Cas9 fragment 1 using the GeneJet Gel extraction Protocol</p> <hr/> <p>PCR PURIFICATION Purification was done on Cre-Cas9 fragment 2 and 4 using PCR purification protocol. Concentration was measured with Nanodrop.</p>	<p>The following was observed in each well</p> <ol style="list-style-type: none"> 1: OK -> gel electrophoresis 2: OK-> PCR purification 3: no band. -> redo with annealing T 60 C. 4: OK -> PCR purification 5: no band -> redo with annealing T 60 C. <hr/> <p>Concentration after gel extraction of Cre-Cas9 fragment 1 was: c(cre-cas9 1)= 0.3 ng/µl.</p> <p>We decided to discard it and redo the PCR for fragment 1.</p> <hr/> <p>Concentration after PCR purification c(cre-cas9 2)=23.6 ng/µl c(cre-cas9 4)=31,8 ng/µl</p> <p>Good concentrations!</p>

Date: 27/6	Performed by: Group 1
What was done	Observations/comments
PCR amplification of Cre-Cas9 constructs	
Preparation and PCR	
<p>PCR mix was done for fragments 1, 3 and 5 using Phusion PCR protocol.</p> <p>The PCR was done using 60 C as annealing temperature and 1 min elongation.</p>	
Gel electrophoresis	
<p>The fragments were loaded on a gel in following order: 1, 3, 5. the gel was run for 35 min at 75 V, 160 mA.</p> <p>Picture was taken.</p>	<p>The following was observed in each well:</p> <p>1: OK -> Gel extraction</p> <p>3: OK -> Gel extraction</p> <p>5: No band → redo</p>
Gel extraction	
<p>First the sample was prepared for gel extraction. 8 µl loading dye was added to 48 µl PCR-product. It was then loaded on the gel.</p> <p>Gel extraction was done on Cre-Cas9 fragment 1 and 3 using the GeneJet Gel extraction Protocol</p> <p>The concentrations were then measured with Nanodrop</p>	<p>Concentrations for Cre-Cas9 fragment 1 and 3.</p> <p>c(Cas9_1)=23.2 ng/µl</p> <p>c(Cas9_3)=11,2 ng/µl</p> <p>Both were stored in the freezer! :)</p>

Date: 27/6	Performed by: Group 1
What was done	Observations/comments
PCR amplification of Cre-gRNA constructs	
Preparation and PCR	
PCR mix was done for fragments 1, 3, 4 and	

5 using Phusion PCR 'protocol.	
The PCR was done using 60 C as annealing temperature and 1 min extension.	
Gel electrophoresis	
<p>PCR . the gel was run for 35 min av 75 V, 160 mA.</p> <p>Picture was taken.</p>	<p>Following was observed for each fragmentl:</p> <p>gRNA 1: wrong -> redo</p> <p>gRNA3: two bands -> gel extraction</p> <p>gRNA 4: wrong ->redo</p> <p>gRNA5: wrong -> redo</p> <p>After gel extraction the concentration för gRNA 3 was too low.</p> <p>c(gRNA_3)=1.6 ng/µl</p>
Date: 27/6	Performed by: Group 1
What was done	Observations/comments
PCR amplification of Cre-gRNA constructs and Cre-Cas9_5 construct	
Preparation and PCR	
<p>PCR mix was done for gRNA fragments and 1, 4 and 5 and also Cas9_5 using Phusion PCR protocol. The PCR was run over night.</p> <p>The PCR was done using 60 C as annealing temperature and 1 min extension.</p>	
Gel electrophoresis	
<p>The day after the gel was run for 40 min av 75 V, 160 mA.</p> <p>Picture was taken.</p>	<p>Following was observed.</p> <p>gRNA 1: weak band at ~1000 bp -> PCR purification</p> <p>gRNA 4: no band -> redo</p> <p>gRNA 5: weak band at ~400 bp -> PCR purification</p> <p>Cas9_5: no band. -> troubleshoot and redo!</p>
PCR purification of gRNA fragments 1 and 5.	
PCR purification was done according to protocol.	Concentrations was:

After the purification the concentrations were low and we decided to redo the PCR again.

c(gRNA1)=3,2 ng/μl
c(gRNA5)=2,0 ng/μl

Date: 27/6

Performed by: Group 1

What was done

Observations/comments

Fusion PCR of CreCas9s construct

PCR

The fusion PCR were done according to the Fusion PCR protocol.

Mixture	size	Total ng
4.5+m		
4+4.5+m		

Mixture	Sample 4 [ng]	Samp 4.5 [ng]	Sample m [ng]
4.5+m			
4+4.5+m			

Mixture	4 [μl]	4.5 [μl]	m [μl]

It was then loaded on a gel.

The gel picture was unclear and smudgy but this is what we got from it:


They seemed to work.	4.5+m seems to have a clear band at 4300 bp. -> PCR purification
	4+4.5+m had one strong band at 4800 bp and one weak at 1500 bp -> gel extraction
PCR purification and gel extraction was done on the Fusion PCR fragments.	c(Cas9_4+4.5+m)=2.4 ng/μl Low, but we could use this and amplify this fragment!!
The concentration for 4+4.5+m after gel extraction was 2.4 ng/μl.	c(Cas9_4.5+m)=16 ng/μl Ok concentration. Could use this and amplify this fragment!
The concentration for 4.5+m after PCR purification was 16 ng/μl.	We will do one last gel electrophoresis to really confirm that the length of each fragment is correct!! This is done on thursday 29/6.

Wednesday 28/6

Date: 28/6	Performed by: Group 1																											
What was done	Observations/comments																											
Fusion PCR of CreCas9s construct																												
PCR																												
<p>The fusion PCR were done according to the Fusion PCR protocol.</p> <p>In the table below is the fragments (Cas9_1-Cas9_4) that will be fused together, and the size of them, and the total amount loaded.</p> <table><tr><td>Mixture</td><td>size [kb]</td><td>Total ng</td></tr><tr><td>1 & 2</td><td>2</td><td>100</td></tr><tr><td>1 & 2 & 3</td><td>3</td><td>150</td></tr><tr><td>3 & 4</td><td>1,5</td><td>95</td></tr></table> <table><tr><td>Mixture</td><td>Sample 1 [ng]</td><td>Samp 2 [ng]</td><td>Samp 3 [ng]</td><td>Samp 4 [ng]</td></tr><tr><td>1 & 2</td><td>50</td><td>50</td><td>-</td><td>-</td></tr><tr><td>1 & 2 & 3</td><td>30</td><td>90</td><td>30</td><td>-</td></tr></table>	Mixture	size [kb]	Total ng	1 & 2	2	100	1 & 2 & 3	3	150	3 & 4	1,5	95	Mixture	Sample 1 [ng]	Samp 2 [ng]	Samp 3 [ng]	Samp 4 [ng]	1 & 2	50	50	-	-	1 & 2 & 3	30	90	30	-	<p>Amanda did a new protocol, which is much easier to follow. If something is unclear talk to her.</p>
Mixture	size [kb]	Total ng																										
1 & 2	2	100																										
1 & 2 & 3	3	150																										
3 & 4	1,5	95																										
Mixture	Sample 1 [ng]	Samp 2 [ng]	Samp 3 [ng]	Samp 4 [ng]																								
1 & 2	50	50	-	-																								
1 & 2 & 3	30	90	30	-																								

3 & 4	-	-	35	35
Mixture	Sample 1 [μl]	Sample 2 [μl]	Sample 3 [μl]	Sample 4 [μl]
1 & 2	2,2	1,95	-	-
1 & 2 & 3	1,3	3,8	2,7	-
3 & 4	-	-	3	1

The gel results from first round fusion PCR:



The fragments 1&2 seems to have fused correctly

Fragments 1&2&3 produced a faint band at the correct position when fused

Fragments 3&4 did not produce a band after the first round of fusion PCR

Second round fusion PCR of Cas9 fragments

All Cas9 fragments (12, 123, 34) were amplified again using their respective end primers. The second round PCR was done according to the fusion PCR protocol Amanda made.

Fragment 12 used the primers 1F and 2R


Fragment 123 used the primers 1F and 3R

Fragment 34 used the primers 3F and 4R

PCR was left overnight

Date: 26/6	Performed by: Group 1
What was done	Observations/comments
PCR amplification of Cre-gRNA constructs	
Preparation and PCR	
<p>PCR mix was done for gRNA fragments and 2 and 4.5 and also Cas9_5 using Phusion PCR protocol. The PCR was run over night.</p> <p>The PCR was done using 55 C as annealing temperature and 1 min extension.</p>	
Gel electrophoresis	
<p>The the gel was run for 40 min av 75 V, 160 mA.</p> <p>Picture was taken.</p> <p>gRNA_2 was okay! :D</p>	<p>Following was observed.</p> <p>gRNA 2: band at ~1000 bp -> PCR purification</p> <p>gRNA 4.5: no band -> redo using PrimeSTAR.</p> <p>If PrimeSTAR does not work, we could try to to a Touchdown PCR.</p>
PCR Purification of gRNA_2	
<p>Purification was done according to protocol.</p> <p>The concentrations was good, the fragment could be used for Fusion!</p>	c(gRNA_2)=15.3 ng/μl

Date: 26/6	Performed by: Group 1
What was done	Observations/comments
PCR amplification of Cre-gRNA constructs	
Preparation and PCR	
<p>PCR mix was done for gRNA fragments1, 3, 4 and 5 and also Cas9_5 using PrimeSTAR PCR protocol.</p>	<p>Apperently when you are using PrimeSTAR Polymerase you should use 55 C annealing temperature. That is the optimal temperature!!!</p>

The PCR was done using 60 C as annealing temperature and 1 min extension.	And also respect the extension time!
Gel electrophoresis	
<p>The fragments were loaded on a gel! 75 V, 40 min, 160 mA.</p> 	<p>Observed:</p> <p>gRNA1 had a very faint band but at the correct length. Because the band was so faint no purification was done. We decided to redo the PCR.</p> <p>gRNA3 had a band at ~900, correct. PCR purification.</p> <p>gRNA4 had a band ~700, correct. PCR purification.</p> <p>gRNA5 had a band at ~300, correct. PCR purification</p>

Thursday 29/6

Date: 29/6	Performed by: Group 1
What was done	Observations/comments
<p>PCR amplification of Cre-gRNA constructs 1,3 and 4.5 and Cre-Cas9 construct 5</p>	
Preparations and PCR	
The PCR mixture was prepared using PrimeSTAR PCR Protocol for gRNA fragments 1, 3 and 4.5 and Cas9 fragment 5.	

The annealing temperature was 55 C and the extension time was 1 min.

Gel electrophoresis

gRNA fragment 1 produced a very weak fragment around the expected 1kb mark. Too weak and unspecific ⇒ redo!

gRNA fragment 3 produced correct bands at 1kb, however the yield was very low compared to the other smaller sized bands ⇒ redo!


gRNA fragment 4.5 produced no band.

Cas9 fragment 5 produced a very strong band at about 500kb which is the expected size! It was PCR purified at a concentration of 141 ng/μl.

The fusion PCR Cas9 fragments 12 and 123 produced no band.

A gel electrophoresis was done with gRNA fragments 1, 3 and 4.5 and Cas9 fragment 5 from the previous PCR.

We also loaded the Cas9 fusion PCR fragments 12 and 123 from the 28/6 to see if the gel loaded correctly.



Date: 29/6

Performed by: Group 1

What was done

Observations/comments


PCR amplification of Cre-gRNA constructs 1,3 and 4.5

Preparations and PCR	
<p>The PCR mixture was prepared using Phusion PCR Protocol for gRNA fragments 1, 3 and 4.5.</p> <p>The annealing temperature was 58 C and the extension time was 1 min.</p>	
Gel electrophoresis	
<p>A gel electrophoresis was done on gRNA fragments 1, 3 and 4.5 @75V for 45mins.</p> <p>No gel image available, however the results are as following:</p> <p>1: no band</p> <p>3: multiple bands, one at correct position ⇒ gel extraction!</p> <p>4.5: no band</p>	
Date: 29/6	Performed by: Group 1
What was done	Observations/comments
Fusion PCR of CreCas9s construct	
Gel electrophoresis	

A gel was loaded with our Cas9 fusion constructs from the 28/6 and 27/6.

These fused constructs are 1&2, 1&2&3, 3&4, 4&4.5&m, 4.5&m

Gel pictures after second round PCR:



Fragments 12 and 123 did not produce any band

Fragments 34 produced a correct band @1.3kb when fused \Rightarrow gel extraction!!

4.5&5 produced 2 bands, however 4&4.5&5 produced a single band \Rightarrow PCR purification was done on that band!

Date: 29/6	Performed by: Alex & Angelo						
What was done	Observations/comments						
<p style="text-align: center;">Gibson assembly for Cas9 with p416 vector</p>							
<p>The Plasmid p416 was cut between H/R regions and gel purified before the gibson assembly.</p> <p>The concentration of the cut plasmid was 12.9 ng/μl.</p> <p>The DNA amounts(μl) for each Cas9 fragment:</p> <table border="1"> <thead> <tr> <th>Fragment</th><th>volume (μl)</th></tr> </thead> <tbody> <tr> <td>1</td><td>1.8</td></tr> <tr> <td>2</td><td>1.6</td></tr> </tbody> </table>	Fragment	volume (μ l)	1	1.8	2	1.6	<p>Due to the difficulty of fusing our fragments using fusion PCR before performing the Gibson assembly, we opted to just performing it without fusion PCR.</p> <p>This involves combining our 7 fragments and our plasmid vector all at once with the Gibson mastermix. It should theoretically work to combine these many fragments, as other researchers have been able to combine even more at a reasonable yield.</p>
Fragment	volume (μ l)						
1	1.8						
2	1.6						

3	3.5	It should be noted that the protocol calls for a 1:1 volume ratio between Gibson mastermix and DNA. We used a 1:2 ratio
4	0.7	
5	0.5	
4.5	1.9	
M	6.2	
Vector	3.9	
Total	20.1	
This was added to an eppendorf tube with 10µl of Gibson mastermix.		

Date: 29/6	Performed by: Group 1
What was done	Observations/comments
PCR amplification of repair fragments O1, O2, R1 and R2	
Preparations and PCR	
<p>The PCR mixture was prepared using Phusion PCR Protocol for O1, O2, R1, R2. These correspond to the following tubes from IDT:</p> <p>upStr2_ptef1_ratI7----->R_1 RatI7_CYCTer_dnSTE2----->R_2 ste3up_olfr1258----->O_1 Olfr1258_ste3dwn----->O_2</p> <p>The annealing temperature was 52 C and the extension time was 1 min.</p> <p>The PCR was left to run overnight.</p>	

Date: 29/6	Performed by: Alex	
What was done	Observations/comments	
Yeast transformation of STE2gRNA and STE3gRNA, and Raphaels Cas9 plasmid with gRNA		
O/N-culture		
2 yeast culture was inoculated in 5 mL YPD along with a negative control (falcon 50 mL tubes). One colony was picked from each strain. The O/N-culture were put in a shake room 30°C overnight.		
Yeast strain	deletion	Mating type
11C	$\Delta HIS, \Delta URA$	a
61A	$\Delta HIS, \Delta URA, \Delta LEU$	α

Friday 30/6

Date: 30/6	Performed by: Group 1
What was done	Observations/comments
Gel Extraction of gRNA fragment 3	
mass of eppendorf tube = 1.00199g mass of tube with gel = 1.4114g \Rightarrow \Rightarrow 410 μ l binding buffer The final concentration was 3,2 ng/ μ l.	Low final concentration, but we are saving this! If we get low concentration a lot of times on this fragment we can evaporate them and dilute, and add to each other.

Date: 30/6	Performed by: Group 1
What was done	Observations/comments
Gel electrophoresis of PCR products O1, O2, R1 and R2	

SHOW RESULTS!	
---------------	--

Date: 30/6	Performed by: Group 1
What was done	Observations/comments
<p>PCR amplification of repair fragments O1, O2, R1 and R2 as well as gRNA fragment 4.5</p>	
<p>A PCR was prepared using the PrimeSTAR protocol. The fragments used were O1, O2, R1 and R2 as well as gRNA fragment 4.5.</p> <p>To the 4.5 fragment, DMSO (1.5 µl) was added, according to protocol.</p> <p>The PCR was run at an annealing temp of 52C and an extension time of 1min.</p>	
<p>Gel electrophoresis on GPCRs and gRNA 4.5</p>	
<p>The gel was added in the order:</p> <p>O1 O2 R1 R2 gRNA_4.5 gRNA_3</p>	<p>gRNA fragment 3 was loaded to double check that we have the right fragment!</p>

Date: 30/6	Performed by: Alex
What was done	Observations/comments
<p>Yeast transformation of STE2gRNA and STE3gRNA, and Raphaels Cas9 plasmid with gRNA</p>	
<p>Transformation</p>	
<p>Since the first transformation didn't seem to work, a new one was made (No colonies after 3</p>	<p>Following OD was measured after 3 h. (Note that the OD measurement was done after the 3 h</p>

Preparations and PCR

PCR was done on gRNA 4.5 using Touchdown Protocol! Phusion Polymerase was used.

Two tubes were prepared, one with DMSO and one without DMSO.

Following was mixed:

30 µl MQ

10 µl Phusion buffer

1.5 MgCl₂

1 µl dNTP

3 µl forward primer

3 µl reverse primer

1 µl template

0.5 Phusion polymerase

In the tube with DMSO:

2 µl DMSO (3 %)

28 µl MQ

PCR machine was set on the touchdown protocol! Annealing T was first set on 65, then it decrease on degree for every cycle. We had 12 cycles.

The protocol is now in our iGEM folder on the PCR computer!

The original protocol is also in a drawer in our lab bench. The annealing T is first

TM+10

and then

TM

:)